

Express Mail No. EL717079998US

PATENT APPLICATION OF
HARRY W. SCHMIDT, DAVID M. NEAL AND ANTHONY
B. MORRIS
ENTITLED
VIAL HANDLING SYSTEM WITH IMPROVED MIXING
MECHANISM

Docket No. T268.12-0048

VIAL HANDLING SYSTEM WITH IMPROVED MIXING MECHANISM

CROSS REFERENCE TO RELATED APPLICATIONS

5 This application claims the priority of
earlier filed co-pending provisional patent
application Nos. 60/188,665, filed March 11, 2000 and
entitled IMPROVED VIAL HANDLING SYSTEM; and
60/188,269 filed March 10, 2000 and entitled WATER
10 AND SOIL AUTOSAMPLER.

BACKGROUND OF THE INVENTION

 The present invention relates to vial
autosamplers of the type used for laboratory
15 automation. More specifically, the present invention
relates to sample mixing within the vial autosampler.

 Vial autosamplers are used to automate
laboratory analyses associated with gas
chromatography, carbon measurement (total carbon and
20 total organic carbon) as well as other types of
analyses. Typically, a vial autosampler has a
storage area adapted to hold a number of vials to be
analyzed. A robotic system generally grasps one of
the vials and transports it from the storage area to
25 an analytical site. Once transported to the
analytical site, the vial contents are sampled and
the appropriate analysis is performed.

 Autosamplers typically use separate
sampling modules for extracting liquid and gas

samples. One example of such an autosampler is described in U.S. Patent No. 5,948,360 to Rao et al. and assigned to Tekmar Company of Cincinnati, Ohio. Liquid sampling typically involves extracting a known
5 quantity of liquid from the vial that is presented to the sampling module of the autosampler, adding a standard to the sample, and transferring the sample to an analytical device. Under certain situations, the specimen must be diluted by a technician by
10 injecting the specimen with a specified volume of methanol or a water-based solution prior to sampling. The extracted sample or methanol extract is then diluted with water prior to analysis by the analytical device.

15 Gas headspace extraction generally involves injecting the specimen with a solvent, such as water, agitating the specimen, and purging the specimen with a gas. Some autosamplers are adapted to perform static headspace extraction while others are adapted
20 to perform dynamic headspace extraction. In static headspace extraction, the specimen is purged from above the specimen and the headspace is removed and transferred to the analytical device. In dynamic headspace extraction, the specimen is purged from
25 underneath the specimen and the head space is removed and then transferred to the analytical instrument. Autosamplers that are capable of performing the above sample extraction include the Precept II and the 7000

HT autosamplers sold by Tekmar-Dohrmann, of Cincinnati, Ohio.

When sample agitation is desired, a stir member, such as a stir bar is generally provided within the vial. The stir member is designed to interact with magnetic fields. Then, a mixing mechanism subjects the stir member to varying magnetic fields. This is typically done by spinning a magnet either beside the vial, or under the vial. In these instances the magnetic field is simply varying, while remaining essentially stationary. This technique is limited in that the strength of magnetic coupling is not optimal. Thus, in samples where agitation may be a challenge, such as soil samples, the mixing mechanism may fail to generate rotation of the stir member, thus reducing the solvent's effects within the sample.

As sample analysis becomes more and more precise, quantifying concentration down into the parts-per-trillion, it becomes increasingly important to provide very effective sample agitation. A more effective mixing mechanism would indeed provide more effective analysis, while possibly reducing cycle times.

25

SUMMARY OF THE INVENTION

A vial autosampler includes a vial cup adapted to contain a vial with a stir member inside. The vial autosampler includes a vial mixing system for agitating

the contents of the vial. The mixing system has an actuator, such as a motor, and a mixing hub that is coupled to the actuator. The mixing hub includes at least one magnetic field source disposed to rotate a
5 magnetic field about the vial cup.

BRIEF DESCRIPTION OF THE DRAWINGS

Fig. 1 is a perspective view of an illustrative automatic vial handling system with
10 which embodiments of the present invention are useful.

Fig. 2 illustrates a mixing system for vial autosamplers.

15 DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENTS

Fig. 1 shows a perspective view of a vial autosampler device 10 in accordance with the invention. The device 10 includes a base unit 12 that includes a vial storage platform area 14, a
20 sampling station 20, and a fluid handling system comprising valves, glasswork, an other fluid handling components. Sampling station 20 receives a vial containing a specimen and extracts a fluid from the vial for further analysis. Finally, device 10
25 includes a central programmable control circuit that accepts user inputs and controls the operation of device 10.

In operation, a vial is selected from vial storage area 14 and transported to an analytical site.

The vial is generally positioned within a vial cup in the sampling module, which lifts the vial such that a resilient septum, generally on top of the vial, is pierced by a stationary needle. In order to obtain
5 the sample, a solvent may be introduced in the vial, and the contents agitated by a mixing mechanism. Once the vial has been sufficiently agitated, the sample is obtained. In some cases, a purge gas is bubbled through the solvent/sample mixture thereby
10 entraining analytes. The purge gas with entrained analytes is recovered and analyzed in accordance with any suitable technique.

Fig. 2 illustrates a mixing system that is useful with vial autosamplers. One particular
15 example of an analysis employing agitation is the analysis of soil. More specifically, EPA (Environmental Protection Agency) test method 5035, requires the addition of a stir member or stir member to a 40 milliliter vial for soil analysis. The stir
20 member is magnetically coupled to an external agitator such that it rotates, or otherwise agitates the content of the vial.

Fig. 2 illustrates vial 532 resting within vial cup 550, and containing stir member 600. Stir
25 member 600 can be any known stir member of suitable geometry, or any other device suitably able to magnetically couple to magnetic fields. Typically, stir members have a magnetic core which is coated with a chemically inert material such as

polytetrafluoroethylene (PTFE) such as Teflon®, available from DuPont. However, as will be described later in the specification the magnetic coupling provided by the apparatus shown in Fig. 2 is generally stronger than that of previous stirring mechanisms. This is because systems that spin a magnet beside the vial have positions where the magnetic poles of the external magnet and the stir member are not aligned. Moreover, in systems where the magnet spins beneath the vial, the poles of the external magnet and the stir member are generally not as close to each other as systems that provide the magnet beside the vial. Due to the enhanced magnetic coupling of embodiments of the present invention, stir members with metal cores can also be used. Such stir members may in fact be easier to produce and have lower costs than magnetic-core stir members. Such cost savings are significant for test methods where stir rods are added to each and every vial, such as tests in accordance with EPA test 5035.

Mixing hub 602 is rotatably mounted beneath vial cup 550 and includes one or more magnets 604. Mixing hub 602 is generally rotated by an actuator such as a motor, or the like, such that magnets 604 rotate about vial cup 550. In the embodiment shown in Fig. 2, two magnets 604 are shown and mixing hub 602 is driven by stepper motor 606 which is coupled to mixing hub 602 through pulleys 608 and belt 610. Those skilled in the art will recognize that a number

of methods can be used to rotate mixing hub 602 including chain drives, gear drives, direct coupling, and other suitable methods. All such methods are expressly contemplated. Magnets 604 are preferably held within mixing hub 602 by end caps 612. By providing one or more magnets which rotate about vial cup 550, the magnetic field is rotated about the vial, unlike prior art approaches where the field source is relatively stationary and simply spins. Even in embodiments where one magnet is used, stirring member 600 will be drawn to that magnet and essentially pulled around the inner diameter of vial 532. Additionally, when two magnets 604 are used, magnets 604 couple to individual magnetic poles of stir member 600 thus creating a stronger magnetic coupling than previous stirring systems. The enhanced magnetic coupling allows stir member 600 to stir the contents of vial 532 in applications when traditional stir systems would prove ineffective.

Another advantage of the stir system shown in Fig. 2, is that it is possible to continuously stir the contents of vial 532 even as elevator 614 causes vial cup 550 to lift vial 532. Such in transit mixing reduces the amount of time that a vial must sit at a mixing station, thus reducing the total analytic time for vial 532 and thereby increasing total system throughput. Additionally, it can be advantageous to cause elevator 614 to move mixing 602

up and down during mixing thereby causing both agitation in the vertical axis as well.

In one embodiment, motor 606 receives energization commands such that the instantaneous rotational speed of mixing hub 602 is varied as a function of the angular position thereof. For example, if the maximum rotational speed is M , and the instantaneous angular position of mixing hub 602 is θ , then the instantaneous rotational speed, i , is related to $M \sin(\theta)$. Although the $\sin()$ function is preferred for profiled mixing, a number of other mathematical functions can be used.

Another feature of the mixing system shown in Fig. 2 is that passageway 616 of hub 602 aligns with hole 618 in the bottom of vial cup 550. Thus, the combination of passageway 616 and hole 618 facilitate needle rinsing when vial 532 is not present since the needle can eject rinse fluid into vial cup 550 which drains out the bottom of passageway 616.

In vial autosamplers with a single analytical site for both solid and liquid containing vials, it is often necessary to change vial temperature during various stages during the analysis. One challenge, however, is that once vial cup 550 is heated, it may be necessary to wait for vial cup 550, and the vial held within, to cool to a selected temperature before further analysis can be performed. Thus, it is advantageous to provide

features which facilitate the cooling of vial cup 550. To that effect, mixing hub 602 is preferably configured to enhance airflow around vial cup 550 when mixing hub 602 rotates. For example, the shapes
5 of magnets 604 can be selected to be relatively thin rectangles thus providing significant surface area to cause airflow around vial cup 550. Further, fins or vanes can also be provided on mixing hub 602 to increase airflow further. In order to enhance
10 cooling even further, thermoelectric devices, such as Peltier devices, can be mounted on vial cup 550 to thermoelectrically cool vial cup 550. Further still, magnets 604 and the additional fins can be configured to provide airflow around the Peltier devices to
15 further enhance cooling of vial cup 550.

Although the present invention has been described with reference to preferred embodiments, workers skilled in the art will recognize that changes may be made in form and detail without
20 departing from the spirit and scope of the invention. For example, although the mixing system is described with respect to permanent magnets rotating about the vial cup, it is expressly contemplated that electromagnets could be located about the vial at
25 different angular position and successively energized in such a way that the magnetic field(s) rotate(s) about the vial cup without actually requiring physical movement.